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REF	CONTENT		Analyseinstrument(er) hvor <b>cobas c</b> pack kan brukes
08056757190	Alkaline Phosphatase acc. to IFCC Gen.2 (1100 analyser)	System-ID 2011 001	<b>cobas c 303, cobas c 503</b>
Nødvendige (men ikke medfølgende) materialer:			
10759350190	Calibrator f.a.s. (12 x 3 mL)	Kode 20401	
05117003190	PreciControl ClinChem Multi 1 (20 x 5 mL)	Kode 20391	
05947626190	PreciControl ClinChem Multi 1 (4 x 5 mL)	Kode 20391	
05117216190	PreciControl ClinChem Multi 2 (20 x 5 mL)	Kode 20392	
05947774190	PreciControl ClinChem Multi 2 (4 x 5 mL)	Kode 20392	
08063494190	Diluent NaCl 9 % (123 mL)	System-ID 2906 001	

## Norsk

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ALP2: ACN 20110

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**Alkaline phosphatase acc. to IFCC Gen.2**

- 2 King EJ, Armstrong AR. Can Med Assoc J 1934;31:376
- 3 Ohmori Y. Über die Phosphomomesterase. Enzymologia 1937;4:217-231.
- 4 Bessey OA, Lowry OH, Brock MJ. A method for the rapid determination of alkaline phosphatase with five cubic millimeters of serum. J Biol Chem 1946;164:321-329.
- 5 Hausamen TU, Helger R, Rick W, et al. Optimal conditions for the determination of serum alkaline phosphatase by a new kinetic method. Clin Chim Acta 1967;15:241-245.
- 6 Schumann G, Klauke R, Canalias F, et al. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37° C. - Part 9. Reference procedure for the measurement of catalytic concentration of alkaline phosphatase. Clin Chem Lab Med 2011 Sep;49 (9):1439-46.
- 7 Guder WG, Narayanan S, Wisser H, et al. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT-Verlag 1996.
- 8 Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-475.
- 9 Breuer J. Report on the Symposium "Drug effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.
- 10 Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.
- 11 Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. Clin Chem Lab Med 2007;45(9):1240-1243.
- 12 Abicht K, El-Samalouti V, Junge W, et al. Multicenter evaluation of new GGT and ALP reagents with new reference standardization and determination of 37 °C reference intervals. Clin Chem Lab Med 2001;39:Special Supplement pp S 346.
- 13 Estey MP, Cohen AH, Colantonio DA, et al. CLSI-based transference of the CALIPER database of pediatric reference intervals from Abbott to Beckman, Ortho, Roche and Siemens Clinical Chemistry Assays: Direct validation using reference samples from the CALIPER cohort. Clin Biochem 2013;46:1197-1219.
- 14 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

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© 2021, Roche Diagnostics
Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim